# **CENTRIFICO**

## Casting Machine Instruction Manual

# Centrifico Casting Machine

## Instruction Manual

### DESCRIPTION

The Kerr Centrifico Casting Machine permits a controlled centrifugal force which produces the soundest castings. Actuated by a heavy-duty steel spring, the machine utilizes a secondary arm pivoted to the straight arm. The secondary arm carries the flask, keeping it in line with the casting force when the arm revolves. As a result, the cast metal travels in a straight line directly into the mold, producing strong, dense castings.

### UNPACKING

The Centrifico Casting Machine carton contains the following: Item 00009 (with accessories) or 00166 (without accessories), Instruction Manual, Accessories for 00009 – 1C, 2C, 3C Flask Cradles (Fig 1); 1C, 2C, 3J, 4C Flasks (Fig 2); Sprue Bases (Fig 3); Large Clay Crucible (Fig 4); Small Clay Crucible (Fig 5); Universal Tongs (Fig 6); Removable Weights (Fig 7); Casting Knob (Fig 8); Crucible Carriage (Fig 9)

### CAUTION:

- Mount the casting machine in a casting well that is at least 14" deep. This prevents molten metal particles from landing on the operator while arm is in motion.
- Attach the well and casting machine base securely to a solid, level, fireproof surface.
- Balance the machine for each casting before flask burnout. This minimizes costly metal spillage and further assures operator safety. After each balancing, securely tighten the top nut located on the machine arm.
- Use dark tinted safety glasses or goggles to protect the eyes from the bright molten metal and torch flame.
- Wear heat-resistant gloves and apron for protection when melting metal and when casting. Always use tongs when handling hot flasks.
- Lift hand and torch up and away from spinning casting machine arm. Always allow arm to spin to a complete stop. Never attempt to stop casting machine manually.
- Casting machine metal components are not designed for direct contact of a torch flame. Repeated exposure will cause metal fatigue and failure during operation.

### INSTALLATION

Mount the casting machine on a solid, level, fireproof work surface surrounded by the Kerr Standard Casting Well, item 15550 (sold separately) (Fig 10).

Allow 14" / 36cm high and 19" / 49cm diameter for Standard Casting Well installation.

Fig 10

Base (Fig 11) should be securely fastened with the Stop Pin (Fig 12) in the Base positioned to the left of the operator's casting position.

### SET UP OF ARM ASSEMBLY AND ACCESSORIES

Place the arm assembly (Fig 16) on the rotor shaft (Fig 13) so that the pin slots in the arm rest on the balancing pin (Fig 14). Tighten the top nut by hand (Fig 15) onto the rotor shaft to hold the arm assembly securely.

The Crucible Carriage (Fig 17) is already attached to the secondary casting arm (Fig 18). This carriage will accept either the Small (1 oz) or Large (2.5 oz) Clay Crucible, depending on the amount of metal to be cast. Slide either crucible in the crucible carriage with the small-end hole against the hole in the backplate (Fig 19).


Fig 16


The unit accepts the 1C, 2C, 3J or 4C Flasks on the secondary arm. To properly center the flasks to receive the cast metal, the 1C Flask Cradle is suspended on the cradle supporting pins (Fig 20) to hold the 1C Flask. The 2C Flask Cradle is suspended on the cradle supporting pins to hold the 2C Flask. The 3C Flask Cradle is suspended on the cradle supporting pins to hold the 3J Flask. The 4C Flask rests on the secondary arm with no cradle.

The unit is supplied with two counterweights: Lock Weight and Permanent Weight (Fig 21). Two Removable Weights are supplied for additional balancing. Counterweights are required on the unit to balance the weight of the invested flask, crucible and the casting metal. When properly balanced, the casting machine arm assembly will revolve smoothly with a minimum of vibration.

### BALANCING THE ARM

It is important that the casting machine be balanced prior to casting. A properly balanced machine minimizes the possibility of metal spilling from the crucible as well as undue machine bearing wear.

*NOTE:* Be sure to balance the arm before burnout as there is no time to work with a heated flask between its removal from the furnace and the casting procedure. Place the invested flask, prior to burnout, in the cradle with the open end (where the wax is exposed) facing the backplate on the crucible carriage. Be sure that the investment has been scraped from the sides of the flask, and it is even and level with the edges. This will assure proper seating of the flask in the machine. Looking through the hole in the clay crucible, the sprue wax should be visible at the flask opening. This allows one to see the direction the metal will travel in the center of the mold during the casting operation. Next slide the crucible, in its carriage, up against the flask and cradle as tightly as possible. Place the measured amount of metal to be used for the casting in the crucible. This end of the casting arm is now ready to balance against the counterweights on the straight arm.

When balancing the arm of the Centrifico Casting Machine, wedge a small object between the straight arm and secondary arm (Fig 22) so as to extend the secondary arm out straight.

Now loosen the top nut (Fig 23) enough so the arm rocks on the rotor shaft. Adjust the counterweights (Fig 24) on the end of the straight arm as necessary to balance the weight of the flask. Turn the weights to move them in or out until the arm balances. Balance is achieved when a slight tap on either end of the arm will cause it to tilt in that direction.

If weights Lock Weight and Permanent Weight are tight, turn them in opposite direction to loosen. The addition of a Removable Weight placed next to the Permanent Weight may be necessary to achieve balance.

After the weights are adjusted, be sure to retighten top nut.

The machine is balanced for this casting only. Since there is no absolute weight factor for the pattern, metal and invested flask the machine should be rebalanced for each individual casting.

### **WINDING THE CASTING MACHINE**

Install the Kerr Cast-R-Knob (Fig 25) to the end of the straight arm. To wind the casting machine for precious metals, hold the Cast-R-Knob and make three complete clockwise turns of the arm. For non-precious metal castings, wind the machine four times.

The heavy duty spring in the machine will require increasingly more force as one winds. The release of this force later causes the centrifugal casting action.

After the arm is wound, pull up on the rubber tipped Stop Pin (Fig 26) in the machine base. Allow the arm to lock in place against the pin. The crucible should now be facing the operator. This is the casting position. When not holding the arm in its wound position the stop pin drops into the machine base.

### **MELTING THE METAL AND POSITIONING THE FLASK**

Wind the machine and set the Stop Pin (Fig 27) while the burnout is taking place. Hold a lighted melting torch pointed down to the metal in the crucible.

Just before the metal enters the melting stages, remove the flask from the furnace, using casting tongs to place the flask in the Flask Cradle (Fig 28) or between the cradle pins for larger flasks. The sprue opening must be toward the crucible. Slide the crucible carriage backplate (Fig 29) flush against the flask, so that when the metal is cast it will leave the crucible and enter the sprue opening in the invested flask.

Caution: Use dark tinted safety glasses or goggles to protect the eyes. Exposure to bright molten metal and torch flame can cause serious eye damage.

### **CASTING**

Continue to expose the metal to the flame until it is completely fluid. Once fluid, continue to hold the flame over the metal. Hold the Cast-R-Knob enough to release the Stop Pin, allowing it to drop into the base of the machine. Raise the torch and release the casting arm, keeping both the torch and your hand away from the spinning casting arm.

Allow the casting arm to continue spinning until it stops itself. If balancing instructions have been followed, it will spin smoothly and with minimum vibration. As the casting arm spins the metal will cool and solidify in the flask. This helps to obtain a dense, sound casting.

Once the arm has stopped, use casting tongs to remove the flask. Consult with your metal supplier for proper setting times prior to water immersion in removal of investment.

### **LIMITED WARRANTY - LIMITATION OF KERR LIABILITY**

Kerr's technical advice, whether verbal or in writing, is designed to assist dentists and laboratories in using Kerr products. The dentist or laboratory assumes all risk and liability for damages arising out of the improper use of Kerr products. In the event of a defect in material or workmanship, Kerr liability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. To take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental, or consequential damages.

EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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### **Kerr Corporation**
Orange, CA 92867 U.S.A.
### **Kerr Australia Pty. Ltd.**
#### **Japan, Inc.**
Manufactured by Kerr Corporation